Casting models

No need to wait, the impression may be poured immediately after disinfection and up to two weeks, provided that the impression is stored at room temperature. The product is compatible with the best plasters on the market and with polyurethane resins.

Cleaning and disinfecting the dispenser

The dispenser can be cleaned, disinfected, and sterilised by immersion in liquid or in a steam autoclave.


Storage

Impressions made with PERFEXIL must be kept at room temperature. PERFEXIL impression materials are guaranteed for 3 years following the date of production provided they are kept at a temperature of between 5° and 27° C.

Galvanization

PERFEXIL impressions may be galvanised with copper or silver.

Times for clinical use

Specifications

	Light	Regular	Light rapide / Super light	Tray Material	Putty Hard / Putty Soft	Putty Soft Rapide / Putty Super Soft Rapide
ISO 4823 ADA 19	Type 3 Light-bodied Consist.	Type 2 Medium-bodied Consist.	Type 3 Light-bodied Consist.	Type 1 Heavy-bodied Consist.	Type 0 Putty Consist.	Type 0 Putty Consist.
Mixing ratio Base :Catalyst	1 :1	1 :1	1 :1	1 :1 5 :1	1 :1	1 :1 5 :1
Mixing time	Auto mixing system	Auto mixing system	Auto mixing system	Auto mixing system	Auto mixing system	30 sec. Auto mixing system
Total working time*	2 min	2 min	1 min 30 sec	1 min 30 sec	2 min	1 min 30 sec
Time in mouth (minimum)	3 min 30 sec	3 min 30 sec	2 min, 30 sec.	2 min, 30 sec.	3 min 30 sec	2 min 30 sec
Setting time *	5 min 30 sec	5 min. 30 sec.	4 min	4 min.	5 min 30 sec	4 min
Strain in compression (min-max)	3-5%	3-5 %	3-5%	3-5%	1-3%	1-3%
Recovery from deformation	> 99,5%	> 99,5%	> 99,5%	> 99,5%	> 99,0%	> 99,0%
Linear dimensional variation after 24h	< 0,20%	< 0,20%	< 0,20%	< 0,20%	< 0,20%	< 0,20%
Hardness «Shore A»	45±2	45±2	45±2	60±2	60±2	60±2

**The times mentioned must be intended from the start of mixing at 23°C. Higher temperatures reduce the times, lower temperatures increase them.*

Presentations

- PERFEXIL PUTTY HARD: (grey) 1 x 280 ml jar of base
- PERFEXIL PUTTY SOFT: (green) 1 x 280 ml jar of catalyst
- PERFEXIL PUTTY SOFT RAPIDE: (blue)
- PERFEXIL PUTTY SUPERSOFT RAPIDE: (turquoise blue)
- PERFEXIL CARTRIDGES TRAY MATERIAL: (kaki)
- PERFEXIL CARTRIDGES LIGHT: (salmon-pink)
- PERFEXIL CARTRIDGES REGULAR: (purple) for dispensing gun
- PERFEXIL CARTRIDGES LIGHT RAPIDE: (orange) 1 box containing 2 x 50 ml cartridges + 12 mixing tips
- PERFEXIL CARTRIDGES SUPER LIGHT: (yellow)
- PERFEXIL TRAY MATERIAL: (kaki) 1 box containing 2 x 380 ml cartridges + 2 locking rings
- PERFEXIL PUTTY SUPER SOFT RAPIDE: (turquoise blue) + 20 dynamic tips
- DYNAMIC TIPS: bag of 50 pieces
- MIXING TIP LOCKERS: bag of 2 pieces
- MIXING TIPS: bag of 48 pieces
- INTRAORAL TIPS: bag of 96 pieces

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Perfexil platinum +

Addition silicone (vinyl polysiloxane) suitable for high precision impressions

Use

PERFEXIL PUTTY HARD - PUTTY SOFT - PUTTY SOFT RAPIDE – PUTTY SUPERSOFT RAPIDE: Vinylpolysiloxane for first impressions with elastic final consistency, recommended for all impression techniques: two-stage (double impression*) or single-stage (simultaneous).

PERFEXIL TRAY MATERIAL: High viscosity vinylpolysiloxane. Recommended for filling impression trays in conjunction with a 'wash' within the PERFEXIL range (CARTRIDGES LIGHT, LIGHT RAPIDE OU SUPER LIGHT) in the single-stage (simultaneous) technique.

PERFEXIL CARTRIDGES REGULAR: Medium viscosity hydrophilic vinylpolysiloxane. Recommended for use in conjunction with PERFEXIL PUTTY or TRAY MATERIAL in the single-stage (simultaneous) technique, or as a monophase material in a special tray.

PERFEXIL CARTRIDGES LIGHT – LIGHT RAPIDE: Low viscosity hydrophilic vinylpolysiloxane. Recommended for use in conjunction with PERFEXIL PUTTY or TRAY MATERIAL in the two-stage (double impression*) or single-stage (simultaneous) technique.

PERFEXIL CARTRIDGES SUPER LIGHT: Very low viscosity hydrophilic vinylpolysiloxane. Recommended for use in conjunction with PERFEXIL PUTTY or TRAY MATERIAL in the two-stage (double impression*) or in the single-stage (simultaneous) technique.

*When using this product in the two-stage (double impression) technique, we recommend that you place a polyethylene spacer sheet on the putty during the first stage of impression taking to improve detail quality in the next stage.

Impression trays: preparation and adhesives

All standard and/or special trays can be used, provided they are rigid and non-deformable.

Proper adhesion of the material to the impression tray is essential to avoid errors, especially during removal from the mouth. To ensure a strong bond, make sure that the adhesive is appropriate for the type and material of the tray.

IMPORTANT: when using special trays, ask the technician if the resin used is self or light curing. For correct use, follow the instructions for each individual adhesive.

Instructions for use

- **Jars:** PUTTY HARD - PUTTY SOFT - PUTTY SOFT RAPIDE – PUTTY SUPERSOFT RAPIDE
 - Take the 2 components out of the tubs using the colour coded dosing spoons.
 - Mix identical proportions of Base and Catalyst; until you obtain a homogeneous mix, with no streaks (about 30 sec.). Then place into the impression tray.

WARNINGS:

- Proportions must be strictly observed in order to have the right setting times: a larger quantity of catalyst does not accelerate the setting time.
- Close tubs immediately after use. Do not invert lids, lower lids and dosing spoons.
- Contact with latex gloves may prevent polysiloxanes from hardening. Avoid direct or indirect contact with latex gloves. Instead, wash your hands well to eliminate all traces of impurities or use vinyl gloves.
- Take care not to stain your clothing as vinylpolysiloxanes are chemically resistant.

- **Cartridges for dispensing gun:** TRAY-MATERIAL - REGULAR - LIGHT - LIGHT RAPIDE - SUPER LIGHT

Assemble the device as shown in the drawings (1-2).

Warning: Before assembling the tip, ensure that the two components (base and catalyst) flow out evenly by applying a light pressure on the dispenser lever and extruding a small amount of material, which must be removed (3). Subsequently, insert the mixing tip and, if required, the intraoral tip in the cartridge. (4)

- **Cartridges for dispensing machine (ratio 5 : 1):** TRAY-MATERIAL – PUTTY SUPER SOFT RAPIDE

- Press in the levers and remove the protective cap.
- Insert the cartridge into the machine as instructed by the manufacturer.
- Ensure an equal length of material is coming out. Clean and remove any excess material.
- Fit and secure the dynamic-static mixing tip.
- Fit the mixing tip locker taking care to engage the clips on the cartridge.
- The 5:1 system cartridge is now ready to use.
- Leave the dynamic-static mixing tip in place after use. To replace it, simply turn the mixing tip locker through 45° and remove it. Then remove the dynamic-static mixing tip.
- Follow the machine manufacturer's instructions to remove and replace the cartridge.

WARNINGS:

- Always dispense equal amounts of base & catalyst before fitting the dynamic-static mixing tip. If you forget to do so, start the machine and wait until the material coming out is uniform in colour, then take away any non-homogenous material before proceeding.
- Only use original 5:1 system cartridges, dynamic-static mixing tips and mixing tip lockers. The use of other mixing tips does not guarantee optimal mixing and can jeopardise the product's characteristics and end results.
- The mixing tip locker cannot be used with cartridges not from 5:1 mixing system.

Warning

In sensitive subjects, vinyl polysiloxanes may cause irritation or other allergic reactions.

Cleaning and disinfection of the impression

After rinsing well under running water, the impression can be immediately disinfected, always follow the instructions of disinfectant manufacturer.